

# TECHNICAL DATA SHEET

## SPORE STRIPS FOR MONITORING DRY HEAT

### PRODUCT CODE

SKU	Population	Packaging	Quantity
BG-106DH	10 <sup>6</sup>	Glassine pouch	100/box

### PRODUCT DESCRIPTION

Spore Strips for monitoring Dry Heat processes consist of:

- A paper strip (6 mm x 38 mm) inoculated with *Bacillus atrophaeus* NRRL B-4418
- Primary packaging in a white glassine pouch



### INTENDED USE

The Spore Strips may be utilized for routine Dry Heat sterilization process efficacy monitoring, sterilizer qualification testing after installation, relocation, malfunctions, major repairs and sterilization process failures for Industrial applications. The Spore Strips are labeled For Industrial Use Only.

### INSTRUCTIONS FOR USE

Place Spore Strips (a minimum of 10 per exposure is recommended) inside representative materials to be sterilized or within the chamber directly. Package or wrap product as usual, if applicable.

Locate the test packages or Spore Strips in areas most difficult to sterilize, as outlined in your specific sterilization validation protocol (usually four corners front, four corners rear, center-center and center-top) or according to standard operating procedure. Run the cycle.

After sterilization or exposure, remove Spore Strips or product from sterilizer.

Aseptically remove the Spore Strip from the primary packaging and transfer to Soybean Casein Digest Broth (SCDB).

Transfer a minimum of one Spore Strip which has not been exposed in a sterilization process as a Positive Control.

**Incubation:** At least one unused tube of culture medium from the same lot should be incubated with the test series as a Negative Control. Place the cultured Spore Strips, the Positive Control and the Negative Control in an incubator set at 30°C to 37°C. Incubate for a minimum of seven days or per a validated reduced incubation period.

**Monitoring:** Examine the Spore Strips daily during incubation. Record observations.

**Interpretation:** Tubes which demonstrate turbidity with an orange pellicle are considered positive for growth of *Bacillus atrophaeus*. Tubes which remain clear and without pellicle formation are considered negative for growth. For unexpected positives, it is recommended that a Gram stain be performed. Gram positive rods are characteristic of the indicator organism.

**Positive Control:** Tube should demonstrate turbidity with an orange pellicle. If the Positive Control does not result in growth, the exposure is considered invalid. Check the conditions during incubation and verify the capability of the medium to support growth.

**Negative Control:** Tube should remain clear. If growth is present, there is a potential for false positives.

# TECHNICAL DATA SHEET

## SPORE STRIPS FOR MONITORING DRY HEAT

### PHYSICAL PROPERTIES

<b>Process</b>	Dry Heat
<b>Strip Dimensions</b>	6 mm x 38 mm
<b>Glassine Dimensions</b>	25 mm x 76 mm
<b>Packaging</b>	100/box

### MONITORING FREQUENCY

For greatest control of sterilized goods, it is recommended that a minimum of ten (10) Spore Strips be included with every load.

### PERFORMANCE CHARACTERISTICS

<b>Population</b>	1.0 to 5.0 x 10 <sup>6</sup> per strip
<b>Purity</b>	No evidence of contamination present in sufficient numbers to adversely affect the finished product.  D-value at 160°C ± 1°C ≥2.5 minutes  The Dry Heat D-value is based on the requirements outlined in ISO 11138-4.
<b>Dry Heat Resistance</b>	Survival – Kill Times Calculated based on the formulas outlined in ISO 11138-1  z-value ≥20°C  The z-value is based on D-values at three temperatures in the range of 150°C to 180°C. Crosstex typically utilizes D-values determined at 150°C, 160°C and 170°C.
<b>Post-Market Criteria</b>	Population: 50% to 300% of certified population D-value: ± 20% of the certified D-value Survival Time: All Spore Strips result in growth at the certified survival time Kill Time: All Spore Strips result in no growth at the certified kill time

### COMPLIANCE

ISO 11138-1 Sterilization of health care products – Biological indicators – Part 1: General requirements


ISO 11138-4 Sterilization of health care products – Biological indicators – Part 4: Biological indicators for dry heat sterilization processes

Crosstex has a validated method for Total Viable Spore Count. Please inquire for the Technical Bulletin entitled *Population Verification of Paper Carrier Biological Indicators* to ensure consistent methodologies are being utilized when performing verification testing.

# TECHNICAL DATA SHEET

## SPORE STRIPS FOR MONITORING DRY HEAT

### STORAGE AND SHELF LIFE



15°C to 30°C



Keep away from sunlight



20% to 70% relative humidity



Keep dry

#### **Shelf Life**

24 months from the date of manufacture



Protect from heat and radioactive sources



Short excursions outside the range of temperature and relative humidity recommended will not impact the performance of the Spore Strips. Do not use damaged Spore Strips. Do not use after the expiration date. Spore Strips contain live cultures and should be handled with care.

### DISPOSAL

Autoclave for not less than 30 minutes at 121°C or per other validated disposal cycle prior to discard.

Crosstex International, Inc.

6789 W. Henrietta Road, Rush, NY 14543 | T 800.860.1888 | E [productorders@crosstex.com](mailto:productorders@crosstex.com) | [www.crosstex.com/industry-industriallife-sciences](http://www.crosstex.com/industry-industriallife-sciences)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.

©2024 Hu-Friedy Mfg. Co., LLC. All rights reserved.

PN 6968 Rev. A - 4/2024