


Capnography Adapter Instructions

The CO₂ Capnography Adapter fits between the Axess™ Nasal Mask and the Axess™ Scavenging Circuit



- 1 Install the Capnography Adapter into the Axess Mask (*adapter can be inserted on either side of the mask*).
- 2 **NOTE:** *The end of the Capnography Extension Tube should be visible inside the Axess Nasal Mask.*
- 3 Insert the Mixed Gas Fitting (*stepped-shaped fitting*) into the Capnography Adapter.
- 4 Twist/attach the sampling line from the Capnography Monitor to the Luer-Lok® Fitting on the Capnography Adapter.
- 5 Discard the CO₂ Capnography Adapter along with the Mask at the end of the procedure.
NOTE: *Single-Use Only* ☒

NOTE: *Capnography Monitoring is the monitoring of carbon dioxide exhaled by the patient. It is sometimes referred to as CO₂ monitoring or ETCO₂ monitoring.*


CROSSTEX®
A CANTEL MEDICAL COMPANY


ACCUTRON™
Nitrous Oxide/Oxygen Sedation

Luer-Lok® registered trademark of Becton Dickinson and Company.

Axess™ is a trademark of Crosstex International, Inc.

Accutron™ is a trademark of Accutron, Inc., a subsidiary of Crosstex International, Inc.

 Accutron (Crosstex Phoenix)
1733 W Parkside Lane
Phoenix, AZ 85027 USA
+1 800.531.2221
www.accutron-inc.com



MT Promedt Consulting GmbH
Altenhofstrasse 80
D-66386 St. Ingbert/Germany
Te. +49 (0) 6894-581020
Fax +49 (0) 6894-581021



0482

CROSSTEX[®]
A CANTEL MEDICAL COMPANY

