



Microorganism Kill Claims

Sanitex Plus® and Sanicloth® Surface Disinfectants

Microorganism	Sanitex Plus® Spray	Sanitex Plus® Wipes	Sani-Cloth® Plus	Super Sani-Cloth®	Sani-Cloth® HB
Methicillin Resistant Staph Aureus (MRSA)	3 mins	5 mins	5 mins	2 mins	10 mins
Hepatitis B Virus	3 mins	2 mins	2 mins	2 mins	10 mins
Hepatitis C Virus	3 mins	2 mins	2 mins	2 mins	-
Acinetubacter Baumunnii	5 mins	-	-	-	-
Adenovirus Type 5	-	-	-	2 mins	-
Aspergillus Niger	5 mins	-	-	-	-
Candida Albicans	5 mins	-	-	2 mins	10 mins
Canine Distemper	-	-	-	-	10 mins
Canine Parvovirus	5 mins	-	-	-	-
Escherichia Coli (E. Coli) 0157:H7	-	5 mins	5 mins	2 mins	10 mins
Herpes Simples Virus, Type I	-	-	-	-	10 mins
Herpes Simples Virus, Type II	30 secs	30 secs	30 secs	2 mins	30 secs
HIV-1	30 secs	2 mins	30 secs	1 min	30 secs
Influenza Strain A2/HK	30 secs	30 secs	30 secs	2 mins	10 mins
Influenza Strain B Virus (Allen Strain)	-	-	-	-	30 secs
Klebsiella Pneumoniae	10 mins	-	-	2 mins	10 mins
Mycobacterium Bovis (TB)	6 mins	5 mins	5 mins	1 min	-
Pseudomonas Aeruginosa ATCC 15442	3 mins	5 mins	5 mins	2 mins	10 mins
Polio I Virus	3 mins	-	-	-	-
Respiratory Syncytial Virus (RSV)	-	1 min	1 min	1 min	-
Rhinovirus	3 mins	-	-	2 mins	-
Salmonella Choleraesuis	3 mins	5 mins	5 mins	2 mins	10 mins
Staphylococcus Aureus ATCC 6538 (Impetigo)	3 mins	5 mins	5 mins	2 mins	10 mins
Trichophyton Mentagrophytes	5 mins	-	-	-	-
Vaccinia Virus	-	-	-	2 mins	10 mins
Vancomycin-Resistant Enterococcus Faecalis (VRE)	5 mins	5 mins	5 mins	2 mins	10 mins

Note: These products do not eliminate the risk of contracting any disease or infection. Follow all instructions and warnings on the use of these products. Must be used in accordance with applicable federal, state or local regulations and guidelines as may be issued from time to time.