



REQUEST FOR QUOTATION CUSTOM STERILIZATION LABELS

To request a quote for your custom label project, please complete the form below and submit the form with a sketch of your preferred label design. Please email the completed form and your label sketch to InfoIndustrialMarkets@crosstex.com.

Please note that minimum order quantities may apply.

PROJECT OPPORTUNITY

Is this project a new or existing label for your company?	New	Existing
Estimated annual usage		
Preferred order frequency or quantity breaks (i.e., monthly, quarterly or semi-annual order quantities?)		
Please provide a target price		

LABEL SPECIFICATIONS

Individual Label Size	Width:	Height:
Is a tear perforation between the labels required?	Yes	No
Does the adhesive need to be permanent or removable/repositionable?	Permanent	Removable/Repositionable
How will the label be applied?	Manually (By Hand)	By Machine
Specify the surface texture to which the label will be applied	Smooth	Textured
Will the label be required to be chemical or abrasion resistant? If known, please specify the chemicals	If by machine, please outline any special considerations we should know:	

LABEL SPECIFICATIONS CONTINUED

Preferred Configuration	<p>Rolls:</p> <ul style="list-style-type: none"> • Core Diameter 1" 3" Other • Copy Position Number (see chart on next page) • Maximum Outside Diameter of each Roll or • Target Number of Labels per Roll <p>Sheets:</p> <ul style="list-style-type: none"> • Size of Individual Sheets • Target Number of Labels per Sheet <p>Individual labels:</p> <ul style="list-style-type: none"> • Packaging preference <p>Other:</p>
-------------------------	---

WHAT PROCESS WILL THE LABEL BE USED IN OR MONITOR EXPOSURE TO?


<u>Process</u>	<u>Preferred Color Scheme of Chemical Indicator</u>
Steam	Blue to Pink Green to Orange/Tan Pink to Dark Brown/Black Blue to Dark Brown/Black White to Dark Brown Black Yellow to Dark Brown/Black
Ethylene Oxide (EO)	Violet to Green Yellow to Brown Yellow to Blue Red to Yellow
Dry Heat	Light to Dark Blue Orange to Brown Pink to Violet (Purple)
Hydrogen peroxide (VH ₂ O ₂)	Violet to Pink Red to Blue Red to Yellow Yellow to Blue Blue to Pink
Formaldehyde	Violet to Green
Gamma/E-Beam	Yellow to Red Green to Violet
Chemical Vapor	Pink to Dark Brown/Black Blue to Black Yellow to Dark Brown/Black White to Dark Brown/Black
Depyrogenation Processes	Pink to Violet/Brown
Other	
Not Applicable	No Chemical Indicator

LABEL SPECIFICATIONS CONTINUED

Label Type	Label Combined with a Chemical Indicator (CI) Custom Indicator, flood coated substrate with no label portion Autoclave or Sterilizable Label with no Chemical Indicator. Label is required to survive the sterilization process but does not indicate exposure via color transition.
For Label Combined with Chemical Indicator, please specify preferred size and position (please send a label sketch with this form)	
Will the label have any printed colors or text? If yes, how many colors? (Does not include chemical indicator, if needed)	Yes. Please provide artwork No
Is any variable information required? (Barcode, sequential numbering, etc.)	Yes. Please provide artwork No
Will the labels be printed using thermal transfer (uses ribbon) or direct thermal (no ribbon used)?	Yes: Thermal Transfer Direct Transfer Laser No

ROLL ORIENTATION (Applies to Rolls of Labels Only)

Please specify required label orientation by selecting the image below. If the required copy position is not listed, please select option #9. If option #9 is selected, a drawing is required which will need to be sent with this form.



9. For Other Roll Orientations Click Here



OTHER REQUIREMENTS

Do you have any other special requirements?

**IMPORTANT:
PLEASE SUBMIT A LABEL SKETCH WHEN SUBMITTING THIS FORM. IT IS IMPORTANT THE DRAWING IS AS DETAILED AS POSSIBLE.**