



DentaPure[®]

Clean Water. Clear Choice.[™]



Your oral health
and protection are
important to us!

That's why we
purify our dental
unit water.





Dental researchers and microbiologists have confirmed that dental unit waterlines can be contaminated with unacceptable levels of bacteria.

For your protection and oral health we are using **DentaPure®** water treatment cartridges on each of our dental units. This advanced water treatment technology was originally developed for NASA and successfully used for over 20 years to ensure that water consumed in space was free of bacteria, viruses and other harmful organisms.

It truly is Rocket Science!





So, my dental office is using Rocket Science to treat the water in the dental unit waterlines, what does this mean for me and how does it work?

DentaPure reduces bacteria in the dental unit waterlines to a level that meets or exceeds the CDC recommendation for safe drinking water.

DentaPure is a cartridge that is attached to the waterlines in the dental unit. As water passes through the cartridge it releases a small amount of the germicide, iodine, into the water. The iodine kills all the harmful bacteria that are in the dental unit waterlines, delivering safe water to you during your dental procedure. **The elemental iodine in DentaPure contains no allergenic proteins** and is a naturally occurring element in our bodies; it is safe for everyone. You may notice a slight medicine taste during your procedure; not everyone can taste it, but if you do, know that it is DentaPure working to treat the water being used during your procedure.

With so much awareness on the environment we thought it important to share with you that DentaPure is completely safe for the environment, there are no caustic chemicals or hazardous waste associated with the use of DentaPure. Made in the USA, the depleted cartridge may be recycled.

Your oral health and protection are important to us! That is why we made the decision to treat our dental unit water with DentaPure, it is effective and most importantly it is safe.

Be sure with DentaPure

CROSSTEX[®]
A CANTEL MEDICAL COMPANY



10 Ranick Road
Hauppauge, NY 11788

Phone: 888.276.7783

Phone: 631.582.6777

email: crosstex@crosstex.com

Learn more at www.dentapure.com

Be sure with **DentaPure**