



PREVENTIVES


twist[®] plus Oscillating Disposable Prophy Angles

twist[®] plus provides highly efficient and ergonomic polishing

twist[®] plus is a unique disposable prophy angle that moves in a 120° reciprocating (back-and-forth) motion instead of a traditional spinning motion. This patented mechanism greatly diminishes frictional heat while reducing the spatter of prophy paste, saliva, blood and other potentially infectious materials. The innovative ergonomic design optimizes comfort as well as visibility during procedures.



- **120° reciprocating motion** allows for continuous-contact polishing with more control and reduced spatter.
- **Small head design** improves visibility and increases maneuverability.
- **Tight cup seal** on the angle's head reduces the chance of snagging soft tissues.
- **Enhanced, fluted cup design** for improved pick-up and interproximal polishing.
- **Tapered and recessed grip design** ensures a secure, comfortable grip.
- **Robust gear mechanism** provides smooth, continuous running.
- **Soft prophy cup flares easily** for efficient subgingival and interproximal cleaning.
- **Tapered brush** is ideal for use with ortho patients for easy access to areas between brackets, wires and bands.



Ordering Information

Description	Quantity	Ref #
Soft cup, gray	100/Bag	TPLUSPASC
Firm cup, white	100/Bag	TPLUSPAFC
Tapered brush	100/Bag	TPLUSPATB

Crosstex[®] and Twist[®] are registered trademarks of Crosstex International, Inc.



Crosstex International
10 Ranick Road
Hauppauge, NY 11788

888-276-7783
crosstex.com