

## TECHNICAL DATA SHEET



**PRODUCT** Crosstex® Articulating Film

**PRODUCT OVERVIEW** Marks high and low occlusion points and provides the most accurate way to mark every point of articulation. Resulting in clear positive marking without smudging on dry or wet fields.

- MATERIALS**
- Plastic Film with wax pigmented coating.
  - Not made with natural rubber latex.

- DETAILS**
- Superior manufacturing process combined with high quality control guarantees clear, positive marking without smudging on wet or dry fields.
  - Color coded packaging and product descriptions make product identification easy.
  - Articulating film available in blue, red, blue/red combo and black/red combo.
  - Available in single sided and double sided straight formats.

**PRODUCT/  
PACKAGING**

Product is available in the following configurations:

Product Code	Description
	All packages contain 25 sheets per book and 12 books per box
TFMBKR	Black/Red Combo, Precisely True, 0.0005"/ 13 µm
TPP	Blue, Thin Precisely True, 0.0005"/ 13 µm
TFMR	Red Precisely True, 0.0005"/ 13 µm
TFMBR	Red/Blue Combo, Precisely True, 0.0005"/ 13 µm

**STORAGE  
REQUIREMENTS**

Store in dry, room temperature location.

**SHELF LIFE**

Single use only. Professional use only.