

## TECHNICAL DATA SHEET



**PRODUCT** Crosstex® Articulating Foil

**PRODUCT OVERVIEW** Marks high and low occlusion points and provides the most accurate way to mark every point of articulation. Resulting in clear positive marking without smudging on dry or wet fields.

**MATERIALS**

- Mirrored Plastic Film with wax pigmented coating.
- Made with Polyethylene Terephthalate (PET).

**DETAILS**

- Superior manufacturing process combined with high quality control guarantees clear, positive marking without smudging on wet or dry fields.
- Color coded packaging and product descriptions make product identification easy.
- Articulating foil available in red, blue and red/blue combination.
- Available in single sided and double sided straight formats.

**PRODUCT/ PACKAGING** Product is available in the following configurations:

<b>Product Code</b>	<b>Description</b>
	All packages contain 25 sheets per book and 12 books per box
TFOBR	Red/Blue Combo, 0.0005"/12 µm
TFOR2	Red Two-Sided, 0.0005"/12 µm
TFOB2	Blue Two-Sided, 0.0005"/12 µm
TFOB	Blue One-Sided, 0.0005"/12 µm

**STORAGE REQUIREMENTS** Store in dry, room temperature location.

**SHELF LIFE** Single use only. Professional use only.