

California Cleaning Product Right to Know Act of 2017 (SB 258)

Ingredient Disclosure Statement

Below is a list of ingredients used in the manufacturing of the product(s) listed below:

Ref. #: JEZ

Description: CrossZyme™ Concentrated Ultrasonic Cleaner

Note: Ingredients are listed in descending order by weight

Ingredient	CAS No.	Purpose
Water	7732-18-5	Solvent
Boric acid NF	10043-35-3	Buffering agent
EDTA solution HS	6381-92-6	Buffering agent
Chembetaine CP 30%	61789-40-0	Surfactant
Amphoterge KJ2	Mixture	Surfactant
Aerosol DT 75	577-11-7	Emulsifier
Biosoft N-91	68439-46-3	Surfactant
Propylene Glycol	57-55-6	Solvent
Protease A 2.5L	9014-01-1	Enzyme
Lemon Fragrance	NA	Fragrance
Yellow D&C 10	NA	Dye
Blue FD&C 1	3844-45-9	Dye
Citric acid	77-92-9	Buffering agent

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.
© 2020 Crosstex International, Inc. Date Issued: 01/13/2020 Document#: IDS-014-R0