

SaniTex Plus™ 2 Cleaner, Sanitizer & Disinfectant

Ready-to-use, EPA registered intermediate level surface cleaner, sanitizer and disinfectant with a low alcohol/dual quaternary ammonium formulation that provides exceptional cleaning and disinfection in 2 minutes.

Broad Spectrum

Intermediate level tuberculocidal, bactericidal, virucidal and fungicidal surface disinfectant effective against 22 microorganisms

Rapid Action

Proven efficacy with a 2-minute contact time for quick turnover

3-in-1 Product

Clean, sanitize and disinfect with one product

Ready-to-Use

Convenient and easy 2-step process with no mixing required

Compatible

Safe on a wide array of clinical contact and housekeeping surfaces



**2-minute
contact time**



**One product
is all you need**



**Fragrance
free**

Ref. #	Description	Quantity
JSSDP2	SaniTex Plus™ 2 Cleaner, Sanitizer & Disinfectant Spray Bottle	24 oz / 710 mL
JSSGP2	SaniTex Plus™ 2 Cleaner, Sanitizer & Disinfectant Gallon Bottle	1-gallon / 3.8 L



When used as directed, suitable for use on surfaces such as:

- Metals (aluminum, brass, copper, stainless steel, chrome)
- Glass
- Sealed granite
- Vinyl
- Rubber (hard, non-porous)
- Laminated surfaces
- Acrylic
- Plastics (PVC, polycarbonate, polypropylene, polystyrene)
- Clear plastics
- Coated surfaces
- Painted surfaces



Perfect for use at the following sites:

- Dental care facilities
- Hospitals
- Medical facilities
- Outpatient surgery centers
- Endoscopy centers
- Dialysis centers
- Long-term care facilities
- Ambulances



Great for use on many high-touch surfaces where cross contamination is a risk:

- Cabinets
- Chairs/Patient chairs
- Countertops
- Face shields
- Eye wear
- Sinks/Faucets
- Doorknobs
- Light switches
- Computer keyboards/Monitors
- Touch-screens
- Telephones
- Dental equipment surfaces
- Patient care fixtures
- Examination room surfaces
- Ambulance interior surfaces
- Operating room surfaces
- Dialysis treatment equipment surfaces
- Endoscope surfaces
- Instrument reprocessing equipment
- and other similar surfaces...

Effective against the following microorganisms in 2 minutes:

Bacteria

- Mycobacterium bovis* - BCG (TB)
- Acinetobacter baumannii* [ATCC 19606]
- Carbapenem Resistant *Klebsiella pneumonia* (CRKP) [ATCC BAA-1705]
- Escherichia coli* with ESBL [ATCC BAA-196]
- Enterobacter aerogenes* [ATCC 13048]
- Listeria monocytogenes* [ATCC 19117]
- Methicillin-Resistant *Staphylococcus aureus* (MRSA) [ATCC 33591]
- Pseudomonas aeruginosa* [ATCC 15442]
- Salmonella enterica* [ATCC 10708]
- Serratia marcescens* [ATCC 14756]
- Staphylococcus aureus* [ATCC 6538]
- Vancomycin Resistant *Enterococcus faecalis* (VRE) [ATCC 51559]

Viruses


- Avian Influenza A Reassortant Virus (H3N2) [ATCC VR-2072]
- Influenza A Virus (Strain A2/Hong Kong) [ATCC VR-544]
- Swine Influenza A Virus (H1N1) [ATCC VR-333]
- Hepatitis B Virus (HBV)
- Hepatitis C Virus (HCV)
- Herpes Simplex Virus Type 2 (HSV-2) [ATCC VR-734]
- Human Immunodeficiency Virus (HIV-1)
- Rhinovirus (RV-14) [ATCC VR-284]
- Rotavirus (Strain WA)

Fungi

- Trichophyton interdigitale* (Ringworm) [ATCC 9533]



Mycobacterium bovis (TB)



Influenza



All products proudly made in the USA except where noted.

All product names are trademarks of Crosstex International Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.