

May 2015

Dear Valued Customer,

Crosstex/SPSmedical is pleased to announce a new version of its RAPICIDE[®] OPA/28 High-Level Disinfectant. The industry's fastest and one of the longest lasting OPA's now comes specially formulated with a purple tint designed to provide easier identification and differentiation from other clear liquids, as well as lower foaming properties. RAPICIDE[®] OPA/28 High-Level Disinfectant continues to be manufactured with the same base materials and to the same specifications as before. The product is FDA cleared and has passed all lab and clinical testing.

Key points to remember during the conversion process:

- This is an automatic upgrade for accounts currently purchasing RAPICIDE OPA/28 High-Level Disinfectant. Customers will begin seeing the new product any time after May 2015.
- Current RAPICIDE OPA/28 High-Level Disinfectant users do not need to change their procedure to use this new formulation.
- It is acceptable to mix the existing and new formulations. The clear-colored RAPICIDE OPA/28 High-Level Disinfectant from current lots can be mixed with, or used to top-off, the purple-colored RAPICIDE OPA/28 High-Level Disinfectant, and vice-versa, without effect on product efficacy.
- Instructions for use, label claims, and key performance indicators including OPA concentration and operating pH are identical to the original formulation.
- The new formulation will come packaged with a purple starburst on the bottle and a label on the top of the case package for easy identification.
- RAPICIDE OPA/28 High-Level Disinfectant is only compatible with RAPICIDE OPA/28 test strips, which should be used prior to each reprocessing cycle to ensure the OPA concentrate meets or exceeds the MRC level.



ORDERING INFORMATION

Product Description	Prior Catalog Number	New Catalog Number	UOM	Change to ordering
RAPICIDE [®] OPA/28 HLD	ML02-0127	Same	4 gallons /case	Automatic Upgrade
RAPICIDE [®] OPA/28 Test Strips	ML02-0137	Same	2 bottles /case 50 strips /bottle	No Change

RAPICIDE OPA/28 High-Level Disinfectant has been validated for compatibility with the full line of MEDIVATORS[®] AERs supporting reusable chemistry and is cleared for manual high-level disinfection.

If you have any questions or need additional information, please do not hesitate to call Crosstex/SPSmedical at (800) 722-1529 or visit www.opa28.com.

RAPICIDE[®] is a registered trademark of Medivators, Inc.

Sincerely,

John Sullivan
Director of Sales
Crosstex/SPSmedical
Customer Support (800) 722-1529
E-mail: info@spsmedical.com
Web site: www.spsmedical.com