

Helping your patients relax



Complete Solutions for Nitrous Oxide/Oxygen Sedation



CROSSTEX[®]
A CANTEL MEDICAL COMPANY

Did you know?

9-15% of patients avoid seeing the dentist due to anxiety and fear?¹

There are many causes of dental anxiety, whether it's a previous negative experience, the feeling of lack of control or fear of possible pain and discomfort during treatment.

Whatever the reason for the anxiety, patients who avoid dental visits can impact the overall profitability of your practice.

By offering nitrous oxide/oxygen sedation, you can help to increase patient flow and provide a better overall experience for your dental patients.

A 2007 survey by the ADA estimated that 70% of dental practices using any form of sedation employed nitrous oxide-oxygen sedation..²

Nitrous oxide is a colorless, odorless to sweet-smelling inorganic gas that was first used in surgical and dental anesthesia in the mid-1800s.

Today, the combination of inhaled nitrous oxide and oxygen, when used appropriately, can be a safe and effective means of managing anxiety in dentistry.



Advantages of nitrous oxide/oxygen sedation includes:

- Minimal side effects
- Solid safety record
- Cooperative and relaxed patients
- Potential for increased efficiency
- Fewer cancellations
- Increased referrals and better reviews
- Possibility to attract new patients

1. US data, Colgate™ Oral Care Center website

2. Center for Scientific Information, ADA Science Institute Topic updated: January 30, 2017

Nitrous Oxide/Oxygen Sedation Solutions to Suit Your Needs



Nitrous oxide/oxygen sedation can provide a valuable additional revenue stream for your practice while offering patients relaxed, comfortable experiences. Accutron is a market leading product portfolio of innovative nitrous oxide/oxygen sedation equipment and accessories, including single-use masks, scavenging systems, manifolds, flowmeter systems and associated components.

Learn more at accutron-inc.com



Centralized



Portable



Chairside

Nasal Masks for Nitrous Oxide/Oxygen Sedation (Single-use)

Axess™ Low Profile Nasal Mask

Description	Quantity	Ref. #
Small		
Unscented, clear	24/box	53037-9
Birthday Bubblegum	24/box	53037-12
Fresh Mint	24/box	53037-16
Medium		
Unscented, clear	24/box	53035-9
Birthday Bubblegum	24/box	53035-12
Fresh Mint	24/box	53035-16
Large		
Unscented, clear	24/box	53034-9
Birthday Bubblegum	24/box	53034-12
Fresh Mint	24/box	53034-16

Not made with natural rubber latex.



- Low profile design provides unencumbered access to the oral cavity
- Fits snugly without tape or uncomfortable tubes protruding into the patient's nostril
- Lightweight, autoclavable scavenging circuit minimizes drag and pull, enhancing patient comfort
- Available in scented and unscented
- Translucent design allows visual monitoring of patient's breathing
- Individually wrapped for greater infection protection
- Lower cost than competing single-use disposable nasal masks
- Single-use minimizes cross-contamination; saves staff time

ClearView™ Classic

Description	Quantity	Ref. #
Pediatric		
Unscented Grey	12/box	33037-9
Outlaw Orange	12/box	33037-10
Sassy Strawberry	12/box	33037-11
Birthday Bubblegum	12/box	33037-12
Groovy Grape	12/box	33037-17
Adult		
Unscented Grey	12/box	33035-9
Outlaw Orange	12/box	33035-10
Birthday Bubblegum	12/box	33035-12
French Vanilla	12/box	33035-14
Fresh Mint	12/box	33035-16
Large Adult		
Unscented Grey	12/box	33034-9
Outlaw Orange	12/box	33034-10
French Vanilla	12/box	33034-14
Fresh Mint	12/box	33034-16

Not made with natural rubber latex.



- Mask-in-mask design minimizes ambient nitrous oxide in operatory
- Clear outer mask allows visual monitoring of patient's breathing
- Soft inner mask creates improved facial seal
- Low profile expands clinical field of view
- Appealing colors and scents engage and relax patients
- Single-use minimizes cross-contamination; saves staff time
- Compatible with other scavenging circuits, including PIP+

Personal Inhaler Plus+™ Multi-use Masks

Description	Quantity	Ref. #
Large		
Unscented Grey	24/pkg.	33015-9
French Vanilla	24/pkg.	33015-14
Perfect Peach	24/pkg.	33015-15
Fresh Mint	24/pkg.	33015-16
Medium		
Unscented Grey	24/pkg.	33016-9
Outlaw Orange	24/pkg.	33016-10
Sassy Strawberry	24/pkg.	33016-11
Birthday Bubblegum	24/pkg.	33016-12
Piña Colada	24/pkg.	33016-13
French Vanilla	24/pkg.	33016-14
Perfect Peach	24/pkg.	33016-15
Fresh Mint	24/pkg.	33016-16
Small		
Unscented Grey	24/pkg.	33017-9
Outlaw Orange	24/pkg.	33017-10
Sassy Strawberry	24/pkg.	33017-11
Birthday Bubblegum	24/pkg.	33017-12
Groovy Grape	24/pkg.	33017-17

Not made with natural rubber latex.



- Sets the standard when it comes to conscious sedation therapy
- Provides patients with reliable and comfortable nitrous oxide delivery
- Single-use minimizes cross-contamination; saves staff time
- Lightly scented and brightly colored nasal masks
- Unscented PIP+™ nasal masks are also available
- Available in small, medium and large

Nasal Mask Scavenging Circuits

Axess™ Scavenging Circuits

Description	Quantity	Ref. #
Axess Scavenging Circuit for In-line Vacuum Controller	Each	53005
Axess Scavenging Circuit Remote Flow System (RFS™) and Accu-Vac™	Each	53006

Not made with natural rubber latex.



- Steam-sterilizable circuit is environmentally friendly and cost-effective
- Slender, lightweight circuit is more comfortable for patients and allows unobstructed access for clinician
- Minimizes drag and pull, so nasal mask stays in place

ClearView™ Scavenging Circuits

Description	Quantity	Ref. #
<i>For Standard Bag Tee</i>		
ClearView™ Scavenging Circuit II (includes Scavenging Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag and complimentary ClearView Single-Use Nasal Hoods)	Each	43003

ClearView™ Scavenging Circuit I (includes Scavenging Circuit with Vacuum Control Valve and complimentary ClearView Single-Use Nasal Hoods)	Each	43005
--	------	-------

For Remote Flow Systems & Accu-Vac™

ClearView™ Scavenging Circuit II (includes Scavenging Circuit Corrugated Tube, Reservoir Bag and complimentary ClearView Single-Use Nasal Hoods)	Each	43004
--	------	-------

ClearView™ Scavenging Circuit I (includes Scavenging Circuit and complimentary ClearView Single-Use Nasal Hoods)	Each	43006
--	------	-------

Not made with natural rubber latex.



- Autoclavable
- X-ray compatible
- Meets NIOSH and ADA flow rate recommendations
- Fits most brands of flowmeters

PIP+™ Scavenging Circuits

Description	Quantity	Ref. #
<i>For Standard Bag Tee</i>		
PIP+™ Scavenging Circuit II (Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag & complimentary PIP+ Masks)	Each	32007

PIP+™ Scavenging Circuit I (Circuit with Vacuum Control Valve, & complimentary PIP+ Masks)	Each	33010
--	------	-------

For Remote Flow Systems & Accu-Vac™

PIP+™ Scavenging Circuit II (Circuit, Corrugated Tube, Reservoir Bag & complimentary PIP+ Masks)	Each	32009
--	------	-------

PIP+™ Scavenging Circuit I (Circuit & complimentary PIP+ Masks)	Each	33008
---	------	-------

Not made with natural rubber latex.



- Autoclavable
- X-ray compatible
- Meets NIOSH and ADA flow rate recommendations
- Flexible system accommodates most brands of flowmeters

ClearView™, Axess™, RFS™, and Accu-Vac™ are trademarks and registered trademarks of Accutron, Inc., a subsidiary of Crosstex International, Inc., a Cantel Medical Company.



Dedicated to making vital contributions to healthcare through high quality solutions, services and education that improve outcomes for dental professionals and patients.