

### Chemical Indicator Labels

*Available for Steam, EO Gas, or Dry Heat sterilization processes.*

#### Description

Chemical indicator labels provide distinct color change when exposed to sterilization cycle.



#### Features

Steam Indicator label is a 3/4" circle label that turns BLUE to BLACK when exposed to Steam sterilization cycles (250°F - 275°F).

EO Gas Indicator label is a 3/4" circle label that turns YELLOW to ORANGE when exposed to EO Gas sterilization cycles (100% or mixture).

Dry Heat Indicator label is a 3/4" circle label that turns from GREEN to BLACK when exposed to dry heat sterilization cycles of 320°F or above.

Storage: Controlled room temperature or refrigeration for best results.

### ORDERING INFORMATION

Cat. No.	Description	Packaging
SIL-1000	Steam Indicator Label	1000 labels/rl
GIL-1000	EO Gas Indicator Label	1000 labels/rl
DIL-1000	Dry Heat Indicator Label	1000 labels/rl