

Isolator™ Plus N95 Surgical Respirator



NIOSH-approved and FDA-cleared flat-fold **N95 Surgical Respirator** with malleable nose and chin pieces for a better fit.

Indicated for use when treating patients with airborne diseases

- Meets guidelines for TB exposure control
- Fluid resistant to penetration by synthetic blood at 160 mmHg
- 99.9% BFE (bacterial filtration efficiency) at 3 microns
- 99.9% PFE (particulate filtration efficiency) at 0.1 micron
- Fluid resistant inner and outer layers
- Chin pocket construction provides increased comfort
- Not made with natural rubber latex

Although the N95 closely resembles a surgical mask, it is a respirator and must be used in accordance with all OSHA regulations regarding respiratory protection.



Ordering Information

Description	Quantity	Ref. #
Isolator™ Plus N95 Surgical Respirator	Box of 28 blue/white striped	GPRN95

For more information call **1-800-Hu-Friedy**, contact your authorized HuFriedyGroup representative, or visit us at [Crosstex.com](https://www.crosstex.com).

Crosstex International, Inc., a Cantel Medical Company
6789 W. Henrietta Road, Rush, NY 14543 USA | 585.359.0130 | 631.582.677 | www.crosstex.com

All product names are trademarks of Crosstex International, Inc., its affiliates or related companies, unless otherwise noted.
© 2020 Crosstex International, Inc. All rights reserved. DLIT01585 Rev A1020

Crosstex is now a proud member of

