



CHEMISTRIES



A product of **Crosstex**[®]

CrossZyme[®] Concentrated Ultrasonic Cleaner



- **Safe, Enzymatic Formulation** – Provides effective cleaning and easy rinsing over a wide range of concentrations and temperatures, is powerfully chelated to improve performance in hard water conditions, is biodegradable, and phosphate and chloride-free components enhance removal of blood and proteinaceous fluids
- **Neutral pH** – Will not damage or harm instruments/endoscopes as it contains rust inhibitors, no harsh alkaline components, is non-corrosive and is safe for use on plastics, stainless steel, glass and aluminum
- **Color/Scent** – Light Green (slightly cloudy)
- **Scent** – Mild Citrus

BY THE NUMBERS

DILUTION

1:20 = General Use

1:10 = Heavy Use

Description	Quantity	Ref. #
<i>Ultrasonic Cleaner</i>	1 Gallon Bottle	JEZ