

# 2021

# Biological Monitoring Quarterly Compliance Calendar

## How often should I perform biological indicator monitoring (spore testing)?

The CDC recommends testing **at least weekly** with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer's directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- Whenever a new type of packaging material or tray is used.
- After training new sterilization personnel.
- After any sterilizer repair or maintenance.
- After any change in the sterilizer loading procedures.

Circle the date each week that BI testing was performed

### JANUARY

S	M	T	W	T	F	S
					1	2
3	4	5	6	7	8	9
10	11	12	13	14	15	16
17	18	19	20	21	22	23
24	25	26	27	28	29	30
31						

### FEBRUARY

S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28						

### MARCH

S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

**FOR STERILITY ASSURANCE:**  
**EVERY LOAD:** Sure-Check™ Sterilization Pouches  
**EVERY DAY:** SteamPlus™ Integrating Indicators  
**EVERY WEEK:** ConFirm™ Biological Monitoring Services

For test status or questions call customer care at: **1-800-819-3336**

[www.CrosstexBMS.com](http://www.CrosstexBMS.com)

All product names are trademarks of Crosstex International, Inc., its affiliates or related companies, unless otherwise noted.

© 2020 Crosstex International, Inc. All rights reserved. DLIT01584 Rev A 0920



2021

# Biological Monitoring

## Quarterly Compliance Calendar

### How often should I perform biological indicator monitoring (spore testing)?

The CDC recommends testing **at least weekly** with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer's directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- Whenever a new type of packaging material or tray is used.
- After training new sterilization personnel.
- After any sterilizer repair or maintenance.
- After any change in the sterilizer loading procedures.

Circle the date each week that BI testing was performed

### APRIL

	S	M	T	W	T	F	S
					1	2	3
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28	29	30		

### MAY

	S	M	T	W	T	F	S
							1
2	3	4	5	6	7	8	
9	10	11	12	13	14	15	
16	17	18	19	20	21	22	
23	24	25	26	27	28	29	
30	31						

### JUNE

	S	M	T	W	T	F	S
			1	2	3	4	5
6	7	8	9	10	11	12	
13	14	15	16	17	18	19	
20	21	22	23	24	25	26	
27	28	29	30				

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

**FOR STERILITY ASSURANCE:**  
**EVERY LOAD:** Sure-Check™ Sterilization Pouches  
**EVERY DAY:** SteamPlus™ Integrating Indicators  
**EVERY WEEK:** ConFirm™ Biological Monitoring Services

For test status or questions call customer care at: **1-800-819-3336**

[www.CrosstexBMS.com](http://www.CrosstexBMS.com)

All product names are trademarks of Crosstex International, Inc., its affiliates or related companies, unless otherwise noted.  
© 2020 Crosstex International, Inc. All rights reserved. DLIT01584 Rev A 0920



# Biological Monitoring Quarterly Compliance Calendar

## How often should I perform biological indicator monitoring (spore testing)?

The CDC recommends testing **at least weekly** with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer's directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- Whenever a new type of packaging material or tray is used.
- After training new sterilization personnel.
- After any sterilizer repair or maintenance.
- After any change in the sterilizer loading procedures.



### JULY

S	M	T	W	T	F	S
				1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31

### AUGUST

S	M	T	W	T	F	S
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31				

### SEPTEMBER

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30		

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

**FOR STERILITY ASSURANCE:**  
**EVERY LOAD:** Sure-Check™ Sterilization Pouches  
**EVERY DAY:** SteamPlus™ Integrating Indicators  
**EVERY WEEK:** ConFirm™ Biological Monitoring Services

For test status or questions call customer care at: **1-800-819-3336** [www.CrosstexBMS.com](http://www.CrosstexBMS.com)



# Biological Monitoring

## Quarterly Compliance Calendar

**How often should I perform biological indicator monitoring (spore testing)?**

The CDC recommends testing **at least weekly** with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer’s directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- **Whenever a new type of packaging material or tray is used.**
- **After training new sterilization personnel.**
- **After any sterilizer repair or maintenance.**
- **After any change in the sterilizer loading procedures.**



**OCTOBER**

S	M	T	W	T	F	S
					1	2
3	4	5	6	7	8	9
10	11	12	13	14	15	16
17	18	19	20	21	22	23
24	25	26	27	28	29	30
31						

**NOVEMBER**

S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30				

**DECEMBER**

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

**FOR STERILITY ASSURANCE:**  
**EVERY LOAD:** Sure-Check™ Sterilization Pouches  
**EVERY DAY:** SteamPlus™ Integrating Indicators  
**EVERY WEEK:** ConFirm™ Biological Monitoring Services

For test status or questions call customer care at: **1-800-819-3336**

[www.CrosstexBMS.com](http://www.CrosstexBMS.com)

All product names are trademarks of Crosstex International, Inc., its affiliates or related companies, unless otherwise noted.

© 2020 Crosstex International, Inc. All rights reserved. DLIT01584 Rev A 0920

