

ClearView™ Adaptor Easy Order Guide

Accutron's ClearView™ Nasal Hoods connect effortlessly to:

- ClearView Scavenging Circuits I (43005 & 43006)
- ClearView Scavenging Circuit II (43003 & 43004)
- Any existing PIP+ Scavenging Circuit
- Any Porter Brown Scavenging Circuit (no adaptor required)
- Matrx Scavenging System (has one smooth tube on each side & cone with two connections) requires adaptor part number 35360



PIP+™ Adaptor Easy Order Guide

Accutron's Personal Inhaler Plus+ Nasal Hoods connect effortlessly to Accutron Scavenging Circuit I (33008 & 33010) and Scavenging Circuit II (32007 & 32009). They also connect easily to any brand of scavenging system. To order an adaptor, please look for these identifying features and note the Accutron order number.

35000	Porter Brown Scavenging System Has 4 tubes – 2 ports on each side, hood within hood bearing the name "Brown."	35350	Matrx Scavenging System Has one smooth, white tube on each side.
35100	Porter Clean Air Scavenging System Has tube inside tube and wings or vacuum sweeps to draw off excess gas.	35400	McKesson Scavenging System Has clear hood and 1 tube on each side.
35325	Matrx Scavenging System Has spiral, white tubing.	35500	Accutron ClearView Scavenging System Has 4 tubes - 2 ports on each side, colored & scented hood within clear outer hood.