

Product Profile

SIL-250, SIS-250, GIL-250, DMS-250, GPS-250R

Chemical Indicators

Available for Steam, EO Gas or Hydrogen Peroxide sterilization processes.


Description

Chemical indicator strips provide distinct color change when exposed to sterilization cycle. Eight inch strips can be used for larger packs, or torn in half at perforation for smaller packs.

Frequency of Use

"An internal CI should be used within each package, tray or rigid container system to be sterilized. The class of CI chosen will depend upon how many critical process variables are to be monitored and how much information is desired about the sterillization process."

- ANSI/AAMI ST79. Comprehensive guide to steam sterilization and sterility assurance in health care facilities.



ORDERING INFORMATION

| Cat. No. | Description | Packaging |
|----------|---------------------------|---------------|
| SIL-250 | 8" Steam Indicator Strip | 250 strips/bx |
| SIS-250 | 4" Steam Indicator Strip | 250 strips/bx |
| GIL-250 | 8" EO Gas Indicator Strip | 250 strips/bx |
| DMS-250 | 8" DUAL™ Monitor Strip | 250 strips/bx |
| GPS-250R | 4" VH2O2 Indicator Strip | 250 strips/bx |

DUAL is a trademark of SPSmedical Supply Corp.