


CHEMICAL INDICATOR TAPE REFERENCE CHART


STEAM – ST



Special Note:

Accept only when indicator turns a dark color when processed. Other color variations may be indicative of suspect sterilization conditions.


VH₂O₂ – HT



Special Note:

Accept only when indicator turns blue when processed. Other color variations may be indicative of suspect sterilization conditions.


EO GAS – GT



Special Note:

Accept only when indicator turns green in EO gas. Other color variations may be indicative of suspect sterilization conditions.

LEAD – LF2 FREE



Special Note:

Accept only when indicator turns a dark color when processed. Other color variations may be indicative of suspect sterilization conditions.