

Chair Sleeves

REF BCCS

INSTRUCTIONS FOR USE

- 1. Hold the chair sleeve by the opening.
- 2. Slide over the chair back until the entire surface area is covered.
- 3. Slide off for removal.



Crosstex International, Inc., a Cantel Medical Company
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | crosstex.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Not made with natural rubber latex.

Made in China
BCCS-IFU-EN / REV B - 06/2021

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.
©2021 Crosstex International, Inc. All rights reserved.