

Technical Data Sheet

Dental Dams



Product	Latex – Powder Free/Low Protein Dental Dams																															
Product Overview	Used to isolate and protect teeth and mouth during a dental procedure																															
Materials	Manufactured from medium modulus natural latex.																															
Details	<ul style="list-style-type: none"> • High tear-strength • 100% Non-chlorinated • Powder-free with less than 5mg/sheet. • Low protein: extractable protein content is below FDA allowable limit of 50µg/sheet (ASTM D5712) 																															
Storage Requirements	Keep dry and away from sunlight																															
Shelf Life	5 years																															
Packaging	<table border="1"> <thead> <tr> <th>Product Code</th> <th>Description</th> </tr> </thead> <tbody> <tr> <td colspan="2">Mint flavored, Green</td> </tr> <tr> <td>19202</td> <td>5 in x 5 in, Thin gauge, 52 Sheets/Box</td> </tr> <tr> <td>19200</td> <td>5 in x 5 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19201</td> <td>5 in x 5 in, Heavy gauge, 52 Sheets/Box</td> </tr> <tr> <td>19402</td> <td>6 in x 6 in, Thin gauge, 36 Sheets/Box</td> </tr> <tr> <td>19400</td> <td>6 in x 6 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19401</td> <td>6 in x 6 in, Heavy gauge, 52 Sheets/Box</td> </tr> <tr> <td colspan="2">Unflavored, Blue</td> </tr> <tr> <td>19102</td> <td>5 in x 5 in, Thin gauge, 52 Sheets/Box</td> </tr> <tr> <td>19100</td> <td>5 in x 5 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19101</td> <td>5 in x 5 in, Heavy gauge, 52 Sheets/Box</td> </tr> <tr> <td>19302</td> <td>6 in x 6 in, Thin gauge, 36 Sheets/Box</td> </tr> <tr> <td>19300</td> <td>6 in x 6 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19301</td> <td>6 in x 6 in, Heavy gauge, 52 Sheets/Box</td> </tr> </tbody> </table>		Product Code	Description	Mint flavored, Green		19202	5 in x 5 in, Thin gauge, 52 Sheets/Box	19200	5 in x 5 in, Medium gauge, 52 Sheets/Box	19201	5 in x 5 in, Heavy gauge, 52 Sheets/Box	19402	6 in x 6 in, Thin gauge, 36 Sheets/Box	19400	6 in x 6 in, Medium gauge, 52 Sheets/Box	19401	6 in x 6 in, Heavy gauge, 52 Sheets/Box	Unflavored, Blue		19102	5 in x 5 in, Thin gauge, 52 Sheets/Box	19100	5 in x 5 in, Medium gauge, 52 Sheets/Box	19101	5 in x 5 in, Heavy gauge, 52 Sheets/Box	19302	6 in x 6 in, Thin gauge, 36 Sheets/Box	19300	6 in x 6 in, Medium gauge, 52 Sheets/Box	19301	6 in x 6 in, Heavy gauge, 52 Sheets/Box
Product Code	Description																															
Mint flavored, Green																																
19202	5 in x 5 in, Thin gauge, 52 Sheets/Box																															
19200	5 in x 5 in, Medium gauge, 52 Sheets/Box																															
19201	5 in x 5 in, Heavy gauge, 52 Sheets/Box																															
19402	6 in x 6 in, Thin gauge, 36 Sheets/Box																															
19400	6 in x 6 in, Medium gauge, 52 Sheets/Box																															
19401	6 in x 6 in, Heavy gauge, 52 Sheets/Box																															
Unflavored, Blue																																
19102	5 in x 5 in, Thin gauge, 52 Sheets/Box																															
19100	5 in x 5 in, Medium gauge, 52 Sheets/Box																															
19101	5 in x 5 in, Heavy gauge, 52 Sheets/Box																															
19302	6 in x 6 in, Thin gauge, 36 Sheets/Box																															
19300	6 in x 6 in, Medium gauge, 52 Sheets/Box																															
19301	6 in x 6 in, Heavy gauge, 52 Sheets/Box																															
Regulatory Information	Product Code: EIE Classification: Class 1, 510(k) Exempt Regulation: 872.6300, Rubber Dam and Accessories Medical Device Listing: D142786 Single-use																															