

Crosstex is now a proud member of HuFriedyGroup



## Technical Data Sheet

### Dental Dams



<b>Product</b>	Latex – Powder Free/Low Protein Dental Dams																															
<b>Product Overview</b>	Used to isolate and protect teeth and mouth during a dental procedure																															
<b>Materials</b>	Manufactured from medium modulus natural latex.																															
<b>Details</b>	<ul style="list-style-type: none"> <li>• High tear-strength</li> <li>• 100% Non-chlorinated</li> <li>• Powder-free with less than 5mg/sheet.</li> <li>• Low protein: extractable protein content is below FDA allowable limit of 50µg/sheet (ASTM D5712)</li> </ul>																															
<b>Storage Requirements</b>	Keep dry and away from sunlight																															
<b>Shelf Life</b>	5 years																															
<b>Packaging</b>	<table border="1"> <thead> <tr> <th>Product Code</th> <th>Description</th> </tr> </thead> <tbody> <tr> <td colspan="2"><b>Mint flavored, Green</b></td> </tr> <tr> <td>19202</td> <td>5 in x 5 in, Thin gauge, 52 Sheets/Box</td> </tr> <tr> <td>19200</td> <td>5 in x 5 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19201</td> <td>5 in x 5 in, Heavy gauge, 52 Sheets/Box</td> </tr> <tr> <td>19402</td> <td>6 in x 6 in, Thin gauge, 36 Sheets/Box</td> </tr> <tr> <td>19400</td> <td>6 in x 6 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19401</td> <td>6 in x 6 in, Heavy gauge, 52 Sheets/Box</td> </tr> <tr> <td colspan="2"><b>Unflavored, Blue</b></td> </tr> <tr> <td>19102</td> <td>5 in x 5 in, Thin gauge, 52 Sheets/Box</td> </tr> <tr> <td>19100</td> <td>5 in x 5 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19101</td> <td>5 in x 5 in, Heavy gauge, 52 Sheets/Box</td> </tr> <tr> <td>19302</td> <td>6 in x 6 in, Thin gauge, 36 Sheets/Box</td> </tr> <tr> <td>19300</td> <td>6 in x 6 in, Medium gauge, 52 Sheets/Box</td> </tr> <tr> <td>19301</td> <td>6 in x 6 in, Heavy gauge, 52 Sheets/Box</td> </tr> </tbody> </table>		Product Code	Description	<b>Mint flavored, Green</b>		19202	5 in x 5 in, Thin gauge, 52 Sheets/Box	19200	5 in x 5 in, Medium gauge, 52 Sheets/Box	19201	5 in x 5 in, Heavy gauge, 52 Sheets/Box	19402	6 in x 6 in, Thin gauge, 36 Sheets/Box	19400	6 in x 6 in, Medium gauge, 52 Sheets/Box	19401	6 in x 6 in, Heavy gauge, 52 Sheets/Box	<b>Unflavored, Blue</b>		19102	5 in x 5 in, Thin gauge, 52 Sheets/Box	19100	5 in x 5 in, Medium gauge, 52 Sheets/Box	19101	5 in x 5 in, Heavy gauge, 52 Sheets/Box	19302	6 in x 6 in, Thin gauge, 36 Sheets/Box	19300	6 in x 6 in, Medium gauge, 52 Sheets/Box	19301	6 in x 6 in, Heavy gauge, 52 Sheets/Box
Product Code	Description																															
<b>Mint flavored, Green</b>																																
19202	5 in x 5 in, Thin gauge, 52 Sheets/Box																															
19200	5 in x 5 in, Medium gauge, 52 Sheets/Box																															
19201	5 in x 5 in, Heavy gauge, 52 Sheets/Box																															
19402	6 in x 6 in, Thin gauge, 36 Sheets/Box																															
19400	6 in x 6 in, Medium gauge, 52 Sheets/Box																															
19401	6 in x 6 in, Heavy gauge, 52 Sheets/Box																															
<b>Unflavored, Blue</b>																																
19102	5 in x 5 in, Thin gauge, 52 Sheets/Box																															
19100	5 in x 5 in, Medium gauge, 52 Sheets/Box																															
19101	5 in x 5 in, Heavy gauge, 52 Sheets/Box																															
19302	6 in x 6 in, Thin gauge, 36 Sheets/Box																															
19300	6 in x 6 in, Medium gauge, 52 Sheets/Box																															
19301	6 in x 6 in, Heavy gauge, 52 Sheets/Box																															
<b>Regulatory Information</b>	Product Code: EIE Classification: Class 1, 510(k) Exempt Regulation: 872.6300, Rubber Dam and Accessories Medical Device Listing: D142786 Single-use																															

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.  
© 2020 Crosstex International, Inc. DLIT01576 Rev A 0820

Crosstex International, Inc.  
 Phone 888-276-7783 | Fax 631-582-1726  
 10 Ranick Road | Hauppauge, NY 11788-4209