

## TECHNICAL DATA SHEET



**PRODUCT** SANI ProZyme® Enzymatic Detergent

**PRODUCT OVERVIEW** SANI ProZyme Enzymatic Detergent is a cleaner designed to remove blood, mucous, tissue and other proteins from delicate instruments. It contains a proteolytic enzyme for effective cleaning of protein-rich soils. SANI ProZyme is not corrosive to instruments and is low foaming. Can be used for both manual and automatic cleaning.

**MATERIALS**

- Glycerin
- Subtilisins (proteolytic enzymes)

**DETAILS**

- Light green liquid concentrate
- Neutral pH – safe and compatible for use on delicate instruments
- Rust Inhibitors – reduces likelihood of rust, tarnish or corrosion
- Rinsing Agents – leaves instruments spot and residue-free
- Low-foaming – compatible with all equipment and pump systems
- Environmentally-friendly – biodegradable and phosphate-free

**PACKAGING** Product is available in the following configurations

Product Code	Description
JED	SANI ProZyme Enzymatic Detergent 1-gallon bottle

**STORAGE REQUIREMENTS** Cool, dry area

**SHELF LIFE** 2 years

Crosstex® and Sani ProZyme® are a registered trademark of Crosstex International, Inc.