

Technical Data Sheet
Crosstex® Articulating Paper



PRODUCT Crosstex® Articulating Paper

PRODUCT OVERVIEW Marks high and low occlusion points and provides the most accurate way to mark every point of articulation. Resulting in clear positive marking without smudging on dry or wet fields.

MATERIALS Paper with wax pigmented coating.

- DETAILS**
- Superior manufacturing process combined with high quality control guarantees clear, positive marking without smudging on wet or dry fields.
 - Color coded packaging and product descriptions make product identification easy.
 - Articulating foil available in red, blue and red/blue combination.
 - Available in single sided and double sided straight formats.

PRODUCT/ PACKAGING Product is available in the following configurations:

Product Code	Description
	12 sheets per book/12 books per box – unless otherwise noted
TPXT	Extra Thin, Blue 0.0015"/38 µm
TPXXT	XX Thin, Blue 0.00125"/32 µm
TPTH	Thick Blue 0.005"/127 µm
TPT	Thin Blue 0.0028"/71 µm
TPBR	Blue Red Combo, 0.0028"/71 µm
TPH	Horseshoe Articulating Paper, Blue Red Combo 0.0035"/89 µm, 12 sheets per book/6 books per box

STORAGE REQUIREMENTS Store in dry, room temperature location.

SHELF LIFE Single use only. Professional use only.