

TECHNICAL DATA SHEET



<u>PRODUCT</u>	Ambi-ance® Ambidextrous Nitrile Gloves, Powder free					
<u>PRODUCT OVERVIEW</u>	Ambidextrous Nitrile Gloves Powder-Free. Gloves are intended for use when conducting medical examinations, diagnostic and therapeutic procedures to protect patient and user from cross contamination/infection.					
<u>MATERIALS</u>	Nitrile (Carboxylate Butadiene Acrylonitrile)					
<u>DETAILS</u>	<ul style="list-style-type: none"> • Maximum comfort and performance with textured fingertips for improved dexterity, enhanced grip and tactile sensitivity. • Soft stretch Nitrile for a comfortable custom fit. • Puncture and tear resistant barrier. • Beaded cuff for extra protection. 					
<u>PRODUCT/ PACKAGING</u>	Product is available in the following configurations					
	Product Code	Description				
	19010	Ambi-ance® Ambidextrous Nitrile Gloves, Extra Small, 200/box				
	19011	Ambi-ance® Ambidextrous Nitrile Gloves, Small, 200/box				
	19012	Ambi-ance® Ambidextrous Nitrile Gloves, Medium, 200/box				
	19013	Ambi-ance® Ambidextrous Nitrile Gloves, Large, 200/box				
	19014	Ambi-ance® Ambidextrous Nitrile Gloves, Extra Large, 200/box				
<u>PRODUCT SPECIFICATIONS</u>	SIZE	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
	Length (mm)	235	235	235	235	235
	Width (mm)	75±5	85±5	95±5	106±5	115±5
	Palm Thickness (mm) *Single Wall	0.06-0.09	0.06-0.09	0.06-0.09	0.06-0.09	0.06-0.09
	Finger (mm)	0.08-0.14	0.08-0.14	0.08-0.14	0.08-0.14	0.08-0.14
	Beading Thickness	≥0.5	≥0.5	≥0.5	≥0.5	≥0.5
	Cuff Thickness (mm)	0.05-0.08	0.05-0.08	0.05-0.08	0.05-0.08	0.05-0.08
<u>STORAGE REQUIREMENTS</u>	Store in a dry cool place; avoid Ozone.					
<u>SHELF LIFE</u>	Three (3) years from the date of manufacture, when above storage conditions are followed.					

Crosstex® and Ambi-ance® are registered trademarks of Crosstex International, Inc.