



Single-Use Disposables

A product of **Crosstex**[®]

SAFE-FLO[®]

FRQUENTLY ASKED QUESTIONS

- SAFE-FLO[®] Saliva Ejectors with unique one-way valves provide a barrier that prevents the backflow described in the CDC's warning on saliva ejectors
- Easy to use – fully assembled saliva ejector with SAFE-FLO valve replaces your current saliva ejector
- The SAFE-FLO valve slightly shortens the exposed length of the saliva ejector, making it less flexible, allowing easy tissue retraction and better access
- Available in white and clear

A product of **Crosstex**[®]



SAFE-FLO[®] Saliva Ejectors



SAFE-FLO SE Valve

- SAFE-FLO SE Valve's unique one-way valve design provides a barrier that prevents backflow
- Can be used with your favorite brand of saliva ejector, allowing you to continue to use your product of choice while protecting your patients against backflow
- The SAFE-FLO SE Valve closes when the suction is stopped or restricted and prevents backflow



SAFE-FLO HV Valve

- For use with HVE tips
- SAFE-FLO HV Valve has the same unique one-way valve design that provides a barrier to prevent backflow
- The SAFE-FLO HV Valve closes when the suction is stopped or restricted and prevents backflow



SAFE-FLO Adapter

- SAFE-FLO Adapter allows you to use the HV suction line with a standard saliva ejector, providing a higher level of suction to the saliva ejector
- Using the HV valve can help reduce the the bioaerosol and control the volume of water in your patient's mouth, while preventing backflow if your patient closes their mouth

Learn More: Crosstex.com/Safe-Flo

Ordering Information:

Description	Quantity	Ref #
SAFE-FLO [®] Saliva Ejector, White	100/Bag	M1000WHBU
SAFE-FLO [®] Saliva Ejector, Clear	100/Bag	M1001CLBU
SAFE-FLO [®] SE Valve	100/Bag	M1010SUSE
SAFE-FLO [®] HV Valve	100/Bag	M1030HVSE
SAFE-FLO [®] Adapter	50/Bag	M1020ADAP

10 Ranick Road
Hauppauge, NY 11788
Tel: (888) 276-7783
Fax: (631) 582-1726
crosstex.com

Not made with natural rubber latex

¹ Center for Disease Control and Prevention. Guidelines for Infection Control in Dental Health Care Settings-2003. MMWR 2003;52 (no. RR-17); [31].



SAFE-FLO[®] is a registered trademark of Crosstex International, Inc., a Cantel Medical Company.

